

North Carolina Division of Air Quality

Quality Assurance Project Plan

2.51 Volatile Organic Compounds (VOCs) by TO-15

Section 4.1

Raleigh Central Office Responsibilities

**Standard Operating Procedure for
Performing a QA Review of a VOC Data Batch**

Version 2015

Submitted by:
North Carolina
Division of Air Quality
1641 Mail Service Center
Raleigh, NC 27699-1641

Approval Sign-Off Sheet

I certify that I have read and approve of the contents of this revision of this SOP with an effective date of 7/31/15.

John Holland 7/23/15
John Holland, UAT Quality Assurance Coordinator Date

Joette Steger 8/11/2015
Joette Steger, PPB Supervisor Date

Corey A. Hughes for Jim Bowyer 8/12/15
Jim Bowyer, LAB Supervisor Date

Jim Bowyer for Donald Redmond 8/12/15
Donald D. Redmond, Jr., Ambient Monitoring Section Chief Date

Table of Contents

2.51.4.1 Standard Operating Procedure for Performing a QA Review of a VOC Data Run

2.51.4.1.1 Purpose

2.51.4.1.2 Equipment Description

2.51.4.1.3 Directions for reviewing the VOC data from the GC/MS laboratory

Addendum: VOC Data QA Checklist

2.51.4.1 Standard Operating Procedure for Performing a QA Review of a VOC Data Run

2.51.4.1.1 Purpose: The purpose of this SOP is to describe the steps needed to review the VOC analysis data from the GC/MS laboratory before it is transferred to the master VOC spreadsheet.

2.51.4.1.2 Equipment Description: PC connected to the network drive group on 'air.ncdenr.net\dfs' – usually mapped to drive letter P.

2.51.4.1.3 Directions for reviewing the VOC data from the GC/MS laboratory:

2.51.4.1.3.1 Open the file folder and check that all the documents on the list below are there. The check list (VOC Data QA Checklist) at the end of this SOP may be used to document this step.

- QA/QC spreadsheet print out and file(s) e-mailed from the VOC Laboratory Analyst
- Leak Check Report
- Sample List Printout
- BFB 524.2 Report
- Calibration Curves Report, may be in a different folder
- Calibration Block Report, may be in a different folder
- VOC Sample Reports in chronological order

2.51.4.1.3.2 Check the dates (the “Acquisition Date” on the VOC Sample Reports) and make sure that all the reports are in chronological order. If there are any missing time periods, check the Sample List Printout and ask the VOC Laboratory Analyst. Each VOC Sample Report (one chromatographic run) should take about 50 minutes and be on 3 pages. The third page may be in a separate group of pages.

2.51.4.1.3.3 Check that the BFB 524.2 Report shows all values “PASS” and that the Leak Check Report shows all inlets end with 0.2 or less for the Entech-Varian system or no fail comments for the Markes-Agilent system. If one or more inlets fail, there may be a subsequent report showing that those leaking have passed.

2.51.4.1.3.4 If there are two QA/QC spreadsheets, combine them into one (see Figure 1). On the QA/QC spreadsheet (see Figure 1) note all the flagged compounds. If a compound of interest like benzene or toluene is flagged, check the details of that flag and void it if possible. For example, if benzene has a flag 3c1 and in checking the calibration data, which may be in a different folder, the r^2 value is 0.9985 and the %RSD is 30.1%, void the 3c1 flag (delete the flag from the QA/QC spreadsheet).

2.51.4.1.3.5 Match the Daily Check Std. Values for the 1st and 2nd runs with the actual VOC Sample Reports and write the date and time of each run on the VOC Data QA Checklist. The first run should be before any samples are run and the second run should be after all samples have been run. These VOC Sample Reports will

- 2.51.4.1.3.7 Review the duplicate sample analyses (Run1 and Run2) in the columns just to the left of the flags columns. Correct the flag 5's if needed.
- 2.51.4.1.3.8 Now review each VOC Sample Report.
- Make sure the system identified at the top is the same system as the calibration data.
 - Make sure there is a correct Sample ID in the "Inj. Notes:" space. This will be matched up to the "VOC Sample Log No." in the master spreadsheet.
 - Does the GC graph look OK, are there strange peaks or lumps? If so, mark it as a bad GC run.
 - Mark out with a single line the data for any compound that has a flag.
 - Circle any value ≥ 0.1 ppb that has not been flagged.
 - Mark with an arrow any value ≥ 1 ppb that has not been flagged.
 - Note how much the CCl₄ value is different from 0.090 ppb.
 - Compare collocated data with primary data and mark any significant differences.
 - Do the concentrations seem reasonable for ambient air at that site; this could indicate that the IDs have been inadvertently switched.
 - Compare the patterns. If two Sample Reports appear very similar, they may have been the same canister but are mislabeled.
- 2.51.4.1.3.9 Send any questions that arose during the review to the VOC Laboratory Analyst. His or her initials should be in the "Operator Name:" space on each VOC Sample Report.
- 2.51.4.1.3.10 When all questions have been resolved, begin recording QA data in the appropriate spreadsheets, for example RepeatsCV2012.xls, DuplicatesCV2012.xls, and CollocatedCV2012.xls. These spreadsheets are located on the P:\ drive in the P:\Toxics\Urban Air Toxics\VOC's Current Year Data\ folder. Record acrolein and benzene precision data. Other compounds may be added.

Type of Data (at least benzene and acrolein)	File to record the data in
Data from collocated and primary samplers run on the same day at the same site (the collocated sample ID will match the primary ID but end in "D")	CollocatedCV201x.xls
Data from two analyses of the same sample in the same batch (the second analysis should have an ID ending with "S")	DuplicatesCV201x.xls
Data from samples that were analyzed in a previous batch (the ID will end in "R")	RepeatsCV201x.xls

- 2.51.4.1.3.11 Once all the data have been recorded in the precision spreadsheets, add lines to the master VOC spreadsheet for these repeat sample analyses and for the duplicate sample analysis (see step 2.51.4.1.3.10 above). A line is added by placing the cursor on the line below the sample ID (VOC Sample Log No.) that

is being repeated, right click, select "Insert," copy the cells from the line with the sample ID that is being repeated up through the Hold Time Flag column and paste them in the inserted line, then change the sample ID in the inserted line by adding "S" or "R" to the end of the sample ID to indicate the expected type of data.

2.51.4.1.3.12 This completes the review. Initial and date the QA/QC Spreadsheet at the bottom and give it to the Laboratory Supervisor for review.

2.51.4.1.3.13 When the Laboratory Supervisor has approved the analytical run and returns the folder, notify the VOC Laboratory Analyst to transfer the data to the P:\ drive for subsequent transfer to the master spreadsheet. (See SOPs for VOCdataFromGCMS_lab and VOCdata2spreadsheetMacroVxx).

VOC Data QA Checklist

For batch run on _____

Before beginning the VOC data review verify that all these documents are present in the data packet:

___ QA/QC spreadsheet print out, see Figure 1 in the “SOP for Performing a QA Review of a VOC Data Run” (and electronic copy of the spreadsheet file e-mailed from the VOC Laboratory Analyst)

___ Leak Check Report

___ “SampleList Printout”

___ BFB 524.2 Report

For a calibration run:

___ Calibration Block Report

___ Calibration Curves Report

For all data batches:

___ VOC Sample Reports matching the “SampleList Printout” and in chronological order

___ Zero air: _____

___ Initial 1 ppb calibration check _____

___ Final 1 ppb calibration check _____

Update Tracking Chart _____

Date of Review: _____ Reviewer: _____