

Preservative Legend (circle above as needed): (A) cool $\leq 6^{\circ}\text{C}$, (B) cool $< 10^{\circ}\text{C}$, (C) 0.008% $\text{Na}_2\text{S}_2\text{O}_3$ [when chlorine is present], (D) H_2SO_4 to pH < 2 , (E) HNO_3 to pH < 2 , (F) HCl to pH < 2 , (G) H_3PO_4 to pH < 2 , (H) 6N NaOH to pH $> 10 < 11$, (I) pH = 9.3-9.7, (J) zinc acetate & NaOH to pH to > 9 , (K) pH 6-9 ascorbic acid [when chlorine is present], (L) EDTA, (M) NaHSO_4 to pH < 2 , (N) 15% EDTA, (O) ammonium sulfate buffer pH = 9.3-9.7, (P) ascorbic acid [when chlorine is present], (R) Lugols, (S) NaHSO_4 to pH < 2 , (V) analyzed within 15 minutes of sample collection, (Z) filtered in field within 15 minutes