

**Preservative Legend (circle above as needed):** (A) cool  $\leq 6^{\circ}\text{C}$ , (B) cool  $< 10^{\circ}\text{C}$ , (C) 0.008%  $\text{Na}_2\text{S}_2\text{O}_3$  [when chlorine is present], (D)  $\text{H}_2\text{SO}_4$  to pH  $< 2$ , (E)  $\text{HNO}_3$  to pH  $< 2$ , (F)  $\text{HCl}$  to pH  $< 2$ , (G)  $\text{H}_3\text{PO}_4$  to pH  $< 2$ , (H) 6N  $\text{NaOH}$  to pH  $> 10$ – $< 11$ , (I)  $(\text{NH}_4)_2\text{SO}_4$  pH=9.3–9.7, (J) zinc acetate &  $\text{NaOH}$  to pH to  $> 9$ , (K) pH 6–9 ascorbic acid [when chlorine is present], (L) leave no headspace, (M) ferrous ammonium sulfate [when chlorine is present], (P) ascorbic acid [when chlorine is present], (R) Lugols, (T) Trizma [when chlorine is present], (Y) analyzed within 15 minutes of sample collection, (Z) filtered in field within 15 minutes using 0.45 $\mu\text{m}$  pore size